


**P 6.4.2**

**Poisson distribution**

P 6.4.2.1 Statistical variations in determining counting rates  
**CASSY-S**



Statistical variations in determining counting rates (P 6.4.2.1)

For each individual particle in a radioactive preparation, it is a matter of coincidence whether it will decay over a given time period  $\Delta t$ . The probability that any particular particle will decay in this time period is extremely low. The number of particles  $n$  which will decay over time  $\Delta t$  thus shows a Poisson distribution around the mean value  $\mu$ . In other words, the probability that  $n$  decays will occur over a given time period  $\Delta t$  is

$$w_{\mu}(n) = \frac{\mu^n}{n!} e^{-\mu}$$


$\mu$  is proportional to the size of the preparation and the time  $\Delta t$ , and inversely proportional to the half-life  $T_{1/2}$  of the radioactive decay.

Using a computer-assisted measuring system, this experiment determines multiple pulse counts  $n$  triggered in a Geiger-Müller counter tube by radioactive radiation over a selectable gate time  $\Delta t$ . After a total of  $N$  counting runs, the frequencies  $h(n)$  are determined at which precisely  $n$  pulses were counted, and displayed as histograms. For comparison, the evaluation program calculates the mean value  $\mu$  and the standard deviation

$$\sigma = \sqrt{\mu}$$

of the measured intensity distribution  $h(n)$  as well as the Poisson distribution  $w_{\mu}(N)$ .

Cat. No.	Description	P 6.4.2.1
559 83	Set of 5 radioactive preparations	1
559 01	End-window counter for $\alpha$ , $\beta$ , $\gamma$ and x-rays	1
524 033	GM-box	1
524 010	Sensor-CASSY	1
524 200	CASSY Lab	1
591 21	Large clip plug	1
590 02	Small clip plug	1
532 16	Connection rod	2
300 11	Saddle base	2
<i>additionally required:</i> PC with Windows 95/NT or higher		1



Measured and calculated Poisson distribution  
 Histogram:  $h(n)$ , curve:  $N \cdot w_{\mu}(n)$